- GRÚA
- LÈVE – PERSONNES
- HOISTS
- PATIENTENLIFTER
- GRUA
- LIFTEN

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MANUAL AND WARRANTY

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INDEX.

INTRODUCTION.
ASSEMBLY.
INITIAL CHECKS.
SAFETY TIPS.
USE AND MAINTENANCE RECOMMENDATIONS.
TECHNICAL SPECIFICATIONS.
CLEANING.
DISPOSAL METHOD.
FAULTS AND REMEDIES.
GUARANTEE.
POSITION OF SLINGS.

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
UK

- A** Leg
- B** Column
- C** Electric actuator
- D** Base opening pedal
- E** Control box + Batteries
- F** Handset
- G** Hoist arm
- H** Hanger
- I** Removable pin
- J** Front wheel without brake
- K** Rear wheel with brake
- L** Base

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INTRODUCTION.

Please read these operating instructions carefully before using your hoist; it contains useful information on safety, handling and maintenance.

Virmedic has wide experience in manufacturing hoists acquired throughout the years. The constant improvements introduced have bolstered our development as a Technical Aids manufacturing company, which is clearly reflected in our wide range of models and applications.

Virmedic hoists comply with the European Directive on Health Products 93/42/CEE and so carry the CE mark. They have been tested at the Biomechanics Institute of Valencia (IBV) and fulfil all the requirements laid down in standard EN ISO 10535: "Hoists for the transfer of disabled persons. Requirements and test methods".

Virmedic has a range of hoists with wheels from 130 to 300 kg, all in electric versions.

These operating instructions have been written in accordance with the current designs of the products. All of them are constantly subject to a process of continuous improvement. Therefore, we reserve the right to modify any technical data that appears in these operating instructions without prior warning. All the data, figures and measurements included in these operating instructions are approximate and must not be considered as strict technical specifications.

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ASSEMBLY.

Every hoist is fully assembled in the factory to check that it works properly and, depending on the version, is dismantled or folded for packaging and transport.

Virmedic hoists are designed to facilitate assembly without the use of any tool.

Remove all the parts of the hoist from the packaging (some of the parts are heavy and you may need someone else to help you).


Put the brake of the rear wheels on and spread the legs to increase the stability of the assembly.

A- FOR DISMOUNTABLE HOISTS.


1. Check that the handle, located on the central tube of the hoist's column is sufficiently loose to allow the base to be fitted.
2. Insert the column into the base.
3. Tighten up the handles.
4. Remove the pin on the hoist arm. Raise the arm until the holes at the end of the actuator is level with that of the securing plates (inserting the actuator end between the securing plate).
5. Re-insert the pin.

B- FOR FOLDABLE HOISTS.


1. Push the column arm assembly up and remove the folding pin from the base of the hoist.
2. Raise the hoist arm-column and replace the folding pin in the base of the hoist.
3. Loosen the handle on the rear of the column and free the hanger.
4. Remove the pin on the hoist arm. Raise the arm until the hole at the end of the actuator is level with that of the securing plates (inserting the actuator end between the two securing plates).

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
INITIAL CHECKS.

Before using your hoist, check that:

- The legs of the hoist open and close correctly. In the case of hoists with mechanical opening these must provide a certain amount of resistance to prevent them from opening or closing accidentally.
- The wheels turn and revolve normally.
- The brakes on the rear wheels work properly.
- The column is fully attached in its usage position.
- The hanger turns and is adequately balanced.
- There is no wear or deformation of the hanger hooks.
- The slings are not damaged or frayed.

ELECTRICAL EQUIPMENT.

1. Before using the hoist for the first time, it is advisable to charge the batteries for at least 5 hours. (B)
2. Check that the connection of the actuator (C) and that of the handset (D) are firmly connected to the control box.
3. Press the up and down buttons (and those opening and closing the legs) to check that the actuator operates correctly.
4. Check that the red emergency stop button (A) is working correctly; if you press it, with the hoist moving up or down, it must stop.

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SAFETY TIPS.

Familiarise yourself with the hoist's devices and how they work before using it for real.

Your hoist is designed for lifting patients, not for any other purpose.

Check that the weight of the patient does not exceed the maximum weight that the hoist can lift or the sling bear.

Check that the end of the actuator is fitted between the plate on the hoist arm and has its pin correctly inserted.

Do not force the controls and devices of the hoist. All of them are easy to use and do not require excessive force to be applied.

Manoeuvre the hoist by pushing on the handles and never push the column, the hoist arm or the patient.

Move the patient with the legs of the hoist in closed position and with the lift arm as low as possible.

The hoists must be handled gently when moving patients, maintaining a speed movement appropriate to the situation.

Move the hoist on flat and straight floors. Their use on surfaces with a slope greater than 5 degrees is not recommended; should you have to pass through a sloping area, we advise you to have a second person to help you.

Never use an electric hoist in a shower.

Do not charge the batteries of a hoist in a bath or shower.

Never leave a patient alone in a hoist.

Virmedic hoists have been designed and manufactured for use with Virmedic slings and accessories; the slings and accessories of other manufacturers have not been tested by Virmedic and so their use cannot be recommended.

Do not use a frayed or worn sling and check that the size and characteristics are suited to the patient.

Adjust the sling in accordance with the sling usage instructions.

Before lifting, always make sure that the belts of the sling are correctly positioned at the end of the hanger; this must be checked when the belts are taut, but before beginning to lift the patient.

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USE AND MAINTENANCE RECOMMENDATIONS.

1. The hoist has swivel wheels; the rear wheels have brakes that stop them from both turning and swivelling. The brake is applied and released with the foot by means of a lever.

We recommend not putting the brake on the wheels when the user is being lifted so that the hoist can move towards the centre of gravity when the patient is lifted from a chair, bed or any other stationary object, whilst at the same time holding it by the handles.

2. The opening of the base, in order to facilitate, for example, the entrance of the wheelchair, is done by pressing with the foot on the appropriate side of the pedal, or electrically by pressing the appropriate button on the handset.

It is advisable to lubricate all the parts of the hoist where movement takes place. This will improve how it operates and prevent possible wear and tear of the parts through daily use.

Prevent dirt accumulating on the axle and the ball bearings of the wheels.

Regularly check that the screws, nuts, etc. are tight and check for slack between the different components of the hoist.

It is not advisable to leave the hoist in places with a high degree of humidity (near the sea, indoor swimming pools, etc.) or in corrosive atmospheres that may affect their components, unless the hoist has been specifically manufactured for this aim.

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ELECTRICAL EQUIPMENT.

The electric actuator has a safety system that disengages the motor if, in the downward movement of the hoist, it comes up against an obstacle that hinders it from descending (parts of the user's body, furniture, ...).

1. The electrical actuator of the hoists for lifting patients of a weight greater than 150 kg has a red emergency lever that allows it, in the event that an electrical fault occurs when the patient is hanging in the air, to be lowered mechanically. A hole on the top of the lever provides access to a screw which, using a flat screwdriver, regulates the speed of descent.

2. The hoists equipped with detachable batteries have an emergency electrical lowering system, situated on the load display screen. This emergency electrical lowering system is operated by inserting the point of a ballpoint pen into that button, which makes the hoist arm descend if a fault has occurred in the handset (do not use a sharp object to press it).

When undertaking the raise-lower movement of the hoist or the open-close movement of the base, keep the appropriate button pressed in, avoiding pressing the buttons repeatedly, as repeated start-ups of the motor uses up more battery, discharging it more quickly.

3. The electrical equipment has an emergency stop button that disconnects the motor's battery automatically halting the movement of the hoist. To disengage the stop button, it must be turned.

When the hoist is NOT in use, press the emergency stop button to prevent the batteries from running down.

Keep the electrical equipment away from water and other liquids to protect it from possible splashes that may affect its internal circuits causing serious damage.

The battery housing has several ventilation holes that ensure the necessary and correct air circulation within this space. This ventilation must not be blocked or covered.

The batteries must be replaced after 4 years and perhaps before, depending on the type of use. Heavy and frequent unloading reduces the life of the battery.

For the optimum life of the battery, recharge it as often as possible. Recharging time is around 5 to 8 hours. If the hoist is used daily, it is a good idea to recharge the battery at night. The control

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box has an automatic device that disconnects it from the charger when it is fully charged. We recommend not undertaking 15 consecutive transfers without charging up the batteries.

Batteries of equipment not habitually used must be recharged at least every 3 months; otherwise they could be damaged as a consequence of their own self-discharge.

To charge, you must check that the connections are correctly attached and that the emergency stop button is not pressed in, otherwise the batteries will not be charged in an appropriate way.

Remember: the battery is the only component of the electrical equipment that requires continuous maintenance. If you do not follow the advice above, the battery may discharge to a level that does not allow it to recover its full charge or may even render it useless.

RECHARGING IN HOISTS WITH FIXED BATTERIES.

The control box has an alarm that sounds when the battery is at approximately 50% of its charge, the moment when it is essential to recharge it immediately. We recommend not waiting to hear the alarm before recharging the battery.

1. To recharge the battery, connect the charger to the connector on the side of the control box and then to the mains supply (220V).

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1. RECHARGING HOISTS WITH REMOVABLE BATTERIES.

The electrical equipment has a screen on which the state of the battery charge is displayed. We recommend that the charge of the battery in use does not fall below 50%, the moment when it is essential to charge it.

2. To recharge the battery you must connect the cable supplied onto the bottom of the control box and then connect it to the main supply (220V). A green light will come on to indicate that the charger is connected to the main supply and a yellow light will remain lit up until the battery is charged.

While the battery is being charged, the hoist cannot be used.

3. The batteries can be recharged using an additional external charger (accessory). To do so, extract the batteries from the hoist by raising the lever on top of them and pull the battery box upwards and outwards.

Place the battery box in the external charger making sure you hear the "click" that confirms they have been fitted properly. Connect the charger to the main supply and then follow the same instructions as above.

TECHNICAL SPECIFICATIONS.

DIMENSIONS & MAXIMUM LOAD.

See the dimensions and weights table on the inside cover of these operating instructions. The dimensions in the table are in centimetres and may vary slightly.

In the height dimensions, the hooks on the ends of the hanger have been taken as the reference point.

The dimensions X, Y and Z have been measured with the hoist legs being open as far as possible.

The Q dimension is the internal width at the position of the hoist's maximum reach (B).

Dimensions S and T of the hoist with sling are a guide and have been taken using a medium back lumbar sling and with the straps in the minimum position.

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ELECTRICAL INFORMATION.

For all its electrical hoist equipment, Virmedic uses the Linak company, a world leader in manufacturing this type of equipment.

The electrical equipment complies with the following standards:

IEC 60601-1 Medical electrical equipment. Part 1: General requirements for basic safety and essential performance.
IEC 60601-1-2 Medical electrical equipment. Part 1-2: General requirements for safety collateral standar: electromagnetic compatibility. Requerements and tests.

The maximum work cycle of the electrical equipment is 10% or 2 minutes of continuous use followed by 18 minutes rest. If this work cycle is exceeded there is a risk that the equipment will overheat and cause a breakdown.

Operational environment of the electrical equipment: from +5° C to +40° C

Actuator 24V continuous current.
Batteries 2 x 12V / 2.9Ah sealed rechargeable batteries of the lead-acid type.

EQUIPO Y PROTECCION IP	Actuador elevación	Actuador apertura	Caja de control	Batería desmontable	Botonera
EQUIPEMENT ET PROTECTION IP	Venin cleration	Venin ouverture	Boîtier de contrôle	Batterie démontable	Télécommande
EQUIPMENT AND IP PROTECTION	Hoist actuator	Open actuator	Control box	Removable battery	Handset
BASICA E-130	LA 31.4 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
BASICA PLUS E-150	LA 3431 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
BASICA XCEL E-175 / XL 175	LA 3431 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
XS 150	LA 31.4 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
HOP E-150	LA 31.4 IP 54	--	CB08 1A IP 30	--	HB71 IP 66
ALEHOP E-175	LA 3431 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
CONFORT E-215	LA 3431 IP 54	--	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66
XXL-300	LA 3431 IP 54	LA 12.1 IP 51	CBJ2 IP 43	BAJ1 IP 65	HB71 IP 66

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ELECTRICAL PROTECTION.

PRODUCT LIFE EXPECTANCY.

All the hoists have been tested in accordance with the standard EN ISO 10535 at the Biomechanics Institute of Valencia (IBV) undertaking a fatigue test over 11,000 raise and lower cycles (hoists), which, depending on the number of hoists per day undertaken, corresponds to an estimated life expectancy of the hoist as indicated in the table below:

		Edad del actuador en años		Age de l'actionneur (années)		Age of actuator in years					
		1	2	3	4	5	6	7	8	9	10
Elevaciones por día	1	365	730	1095	1460	1825	2190	2555	2920	3285	3650
	2	730	1460	2190	2920	3650	4380	5110	5840	6570	7300
	3	1095	2190	3285	4380	5475	6570	7665	8760	9855	10950
	4	1460	2920	4380	5840	7300	8760	10220	11680	13140	14600
	5	1825	3650	5475	7300	9125	10950	12775	14600	16425	18250
	6	2190	4380	6570	8760	10950	13140	15330	17520	19710	21900
	7	2555	5110	7665	10220	12775	15330	17885	20440	22995	25550
	8	2920	5840	8760	11680	14600	17520	20440	23360	26280	29200
	9	3285	6570	9855	13140	16425	19710	22995	26280	29565	32850
	10	3650	7300	10950	14600	18250	21900	25550	29200	32850	36500
Élévations Par jour	11	4015	8030	12045	16060	20075	24090	28105	32120	36135	40150
	12	4380	8760	13140	17520	21900	26280	30660	35040	39420	43800
	13	4745	9490	14235	18980	23725	28470	33215	37960	42705	47450
	14	5110	10220	15330	20440	25550	30660	35770	40880	45990	51100
	15	5475	10950	16425	21900	27375	32850	38325	43800	49275	54750
	16	5840	11680	17520	23360	29200	35040	40880	46720	52560	58400
	17	6205	12410	18615	24820	31025	37230	43435	49640	55845	62050
	18	6570	13140	19710	26280	32850	39420	45990	52560	59130	65700
	19	6935	13870	20805	27740	34675	41610	48545	55480	62415	69350
	20	7300	14600	21900	29200	36500	43800	51100	58400	65700	73000
Hoists Per day	21	7665	15330	22995	30660	38325	45990	53655	61320	68985	76650
	22	8030	16060	24090	32120	40150	48180	56210	64240	72270	80300
	23	8395	16790	25185	33580	41975	50370	58765	67160	75555	83950
	24	8760	17520	26280	35040	43800	52560	61320	70080	78840	87600
	25	9125	18250	27375	36500	45625	54750	63875	73000	82125	91250

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Continue with normal use
Consider changing the actuator
The actuator must be changed.

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The above check list is a guide for LINAK actuators used in conditions of maximum load in accordance with LINAK specifications and in line with the standard EN10535 Patient Hoists.

If the pattern of use, load and run differs from LINAK specifications, this table must be adjusted to suit the application.

The LINAK specifications mentioned can be found on the actuator label, or if it is not there, they can be consulted in the corresponding data sheet at www.linak.com

CLEANING.

The hoist must be cleaned using a damp cloth with a product used for cleaning furniture or similar that does not contain solvents that could damage the coating of the hoist.

Do not allow the electrical elements of the hoist to come into contact with water.

DISPOSAL METHOD.

The hoist is manufactured with steel profiles with an epoxy coating. Therefore, to dispose of it you must follow recycling regulations or local recommendations on disposal.

Rechargeable batteries containing lead (Pb) must be recycled so they should be disposed of in containers destined for recycling cars batteries or at the appropriate waste management sites.

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FAULTS AND REMEDIES.

SYMPTOM

The lift arm, with the hoist unloaded, sometimes does not lower.

CAUSE

· The lift arm needs a minimum weight in order to descend.

SOLUTION

· Lightly push the lift arm downwards at the same time as pressing the down button.

SYMPTOM

The moving parts of the hoist do not function smoothly.

CAUSE

· Lack of lubrication.

SOLUTION

· Lubricate moving parts.

SYMPTOM

The actuator does not work but a "click" can be heard inside the control box when pressing the handset.

CAUSE

- 1 · Actuator plug badly connected to the control box.
- 2 · Flat batteries.
- 3 · Damaged actuator cable.
- 4 · Actuator or control box broken down.

SOLUTION

- 1 · Connect the actuator correctly.
- 2 · Charge batteries.
- 3 · Send actuator to be repaired.
- 4 · Send equipment to be serviced.

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FAULTS AND REMEDIES.

SYMPTOM

The actuator does not work and there is no "click" sound in the control box when pressing the handset.

CAUSE

- 1 · Red emergency button pressed in.
- 2 · Equipment with removable battery that has not made a good contact.
- 3 · Bad connection of the electrical components.
- 4 · Flat batteries.
- 5 · Handset cable damaged.
- 6 · Damaged equipment.

SOLUTION

- 1 · Turn it to release it.
- 2 · Position the battery correctly.
- 3 · Check the connections.
- 4 · Charge the batteries
- 5 · Replace the handset.
- 6 · Send equipment to be serviced.

SYMPTOM

The actuator stops.

CAUSE

- 1 · Flat batteries (acoustic warning sounds on pressing the handset)
- 2 · The weight that the hoist can lift has been exceeded.

SOLUTION

- 1 · Charge batteries.
- 2 · Reduce the weight.

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GUARANTEE.

All hoists have been rigorously inspected; on their labels there is a lot and series number that lists that hoist with the materials, equipment, and components employed and the work processes undertaken for their manufacturing assembly and inspection.

The Virmedic hoists have a three-year warranty for the structure, 5 years for the electrical devices and 1 year for the battery, under normal usage conditions, and respecting the recommendation on use and maintenance and not changing or manipulating any of their components.

Any manipulation or alteration, especially of the electrical equipment renders this guarantee null and void.

In the event of defect or malfunctioning of the hoist or any of its components, **contact the distributor where you purchased it**, presenting it together with the invoice or purchase receipt, so that it may be inspected so that the manufacturer can decide whether to repair it or replace it.

Virmedic accepts no responsibility for imperfections caused by bad usage or as a consequence of not having followed the recommendations given in these operating instructions.

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POSITION OF SLINGS.

1. BACK LUMBAR SLING AND THIGHS SEPARATED.

· Let the back part of the sling fall vertically between the backrest and the patient's back until the material touches the seat.

· The two front points of the sling go under the thighs from outside to inside going round them.

· Fit each of the four straps to the ends of the hanger. In the case of the thigh straps, they may be fitted by crossing them over; that of the right thigh to the left side of the hanger and vice versa.

· It is possible to use it as a legs-together sling if the point that passes under the left thigh also passes under the right thigh and the strap is hung onto the right side of the hanger and vice versa.

2. SEAT-BACKREST LEGS TOGETHER SLING

· To fit it, let the widest part of the sling fall vertically between the backrest of the chair and user's back.

· The user must be lifted and the sling be moved under him/her until it reaches the back of the knee joint.

· On each side of the user there will be two straps that have to be fitted to each end of the hoist's hanger.

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3. SLING TYPE HANGER WITH LEGS SEPARATED.

- To fit it, let the sling fall directly between the chair and the user's back.
- The two front points of the sling are passed under the thighs, from outside to inside going round them.
- Fit each of the four straps to the ends of the hanger. In the case of the thigh straps, these can be fitted by crossing them over; that of the right thigh to the left side of the hanger and vice versa.
- It is possible to use it as a legs-together sling if the point that is passed under the left thigh also passes under the right thigh and the strap is fitted to the right side of the hanger and vice versa.

For all sling types begin to lift the user until the four straps are taut. Choosing the appropriate coloured buckle of each strap, check that all of them are correctly fitted and make sure that the user is in a slightly inclined position and is comfortable and stable before undertaking any movement or transfer.

The handle that some slings have is used to ensure the correct posture of the patient, as well as allowing you to direct the patient's position on the element destined for transfer.

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